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GB User Instructions

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boso

**medicus
s m a r t**

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Your **medicus smart** unit is a state of the art blood pressure computer. It works on the **oscillometric measuring principle**. The oscillations created by the pulse and transmitted via the cuff are stored and evaluated by a micro processor. Therefore there is no need for a microphone the positioning of which is vital if one is to get reliable readings.

The boso-medicus smart conforms to the requirements of the European Directives based on the Law governing Medical Products (CE mark), and to the European Standard EN 1060, Part 1, on „Non-Invasive Blood Pressure Instruments – General Requirements“ and Part 3 „Additional Requirements for Electronic Blood Pressure Systems.“

Regular checks on the accuracy of the measurement system should be carried out in accordance with the valid directives of the individual countries, governing this issue.

The instructions for re-calibration checks are on page 35 of this booklet.

Important Points to observe

- Only use bosso cuffs CA01 and CA02 with the bosso-medicus smart blood pressure monitor.
- Do not expose to extreme temperatures and avoid direct sun light.
- Avoid knocking the unit over or dropping it, and protect it from dirt and moisture.
- Do not place anything on top of the unit and avoid damaging the cuff with sharp instruments such as pins, needles, scissors etc.
- If the unit is not to be used over a period of time the battery should be removed to avoid possible leakage damaging the system.

General Information on Blood Pressure

Your heart, like a pump that works continuously distributing the blood under pressure to every blood vessel throughout your body, ensures that vital, oxygenated blood is delivered to all organs.

Blood pressure is determined by the strength of the heartbeat, the elasticity and the diameter of the blood vessels, as well as various other factors.

The pressure created by the contraction (systolic pressure) of the heart muscles is the higher, and that at the relaxation (diastolic pressure) is the lower blood pressure.

Both these values are necessary for the correct evaluation and a subsequent correct medical diagnosis. Therefore, both the systolic and diastolic pressures are measured when the blood pressure is taken. The measurement unit is "mmHg" = millimeters mercury column.

Blood pressure is never constant, but always changing. In the morning it is lower than in the evening,

and it is lowest when we are asleep. It also changes with every state of physical and mental exertion.

Eating, smoking, fear or stress have an influence on our blood pressure. The seasons too have their effect, with the pressure being slightly lower during summer than during the winter.

In the short term therefore higher blood pressure readings are no cause for concern. If the blood pressure however stays constantly above the average values one speaks of "high blood pressure". Your GP calls this hypertension.

High blood pressure may not immediately be evidenced by any feeling of being unwell. But unobserved and untreated, it can lead to serious health problems. Therefore it is important to have your blood pressure checked by your GP on a regular basis.

With the boso-medicus smart you can now do your own regular check-ups and can considerably contribute to the safeguarding of your health and to the recognition of any disorders at an early stage.

Your blood pressure should be taken twice a day, and if possible, always at the same time: in the morning when you get up and in the evening when the body had a chance to relax again.

Occasional diverging readouts should not cause any undue concern. It is important however to take readings over a number of days and take these into account for a full evaluation.

Blood Pressure (BP) Values

To determine correct blood pressure, the following two values have to be ascertained:

- The **Systolic or upper blood pressure**. This is created when the heart muscle contracts and blood is being forced into the blood vessels.
- The **Diastolic or lower blood pressure**. This is created when the heart muscle extends and once again fills with blood.

Measurement values are expressed in **mmHg (mm mercury column)**.

The World Health Organization (WHO) has set the following guideline for the assessing of blood pressure values:

	Systolic	Diastolic
Too high	over 140 mmHg	over 90 mmHg
Normal – borderline	130 to 139 mmHg	85 to 89 mmHg
Normal	120 to 129 mmHg	80 to 84 mmHg
Optimum	up to 119 mmHg	up to 79 mmHg

Whether medical treatment is necessary depends not only on the blood pressure but also on the risk profile of the individual patient. Please contact your GP if one of the values (SYS/DIA) constantly exceeds the limit (> 140 /> 90).

Important Points to observe when taking your own blood pressure

- Always take your blood pressure in a resting position. Two measurements per day are recommended – in the morning on rising and in the evening when the body is relaxed after the days work.
- Your blood pressure should be taken sitting relaxed and comfortably in an arm chair or at a table where you can support your arm in a slightly angled position.
- The measurement is taken on the free upper arm. Care has to be taken that there is no constriction from pushed-back sleeves (if necessary remove garment).
- The arm has to be kept perfectly still while the blood pressure is taken.
- Between successive measurements an interval of at least 2 minutes is to be allowed.
- The measurement should not last longer than two minutes.
- It is quite normal to experience variations in blood pressure values. Even with repeated tests there can be marked differences. Single or repeated measurements are therefore no indication of the true blood pressure. A reliable evaluation is only possible when you take regular measurements under the same conditions and enter the values in the enclosed test card.
- Taking your own blood pressure does not mean Therapy. Never change any given doses of medication prescribed by your general practitioner.
- Persons suffering any kind of heart defect or those with pace makers are advised to consult their GP, as in these instances altered or incorrect blood pressure readings can occur.
- Using your blood pressure monitor in the vicinity of strong electromagnetic fields (ie Ultra violet units, Mobile telephones) can also cause faulty readings (page 32-34).

Setting-Up Your Unit

Insert Battery

The battery compartment is on the underside of the unit.

- Insert **1 battery (R6 or LR 6)**, 1,5 Volt according to the illustration in the compartment.

Important: If the battery are not inserted in correct polarity it can lead to a build-up of heat and leaking of the batteries.

- Only use high quality leakproof batteries of the given specification (see "Technical Data").
- If the unit is not to be used over a prolonged period of time remove battery.

- The unit shows the following function symbols:

battery fully loaded

battery partially loaded

replace battery soon

flashing! No further measurements possible; replace battery.

- **Care for your environment!**


Used batteries and recharging units do not belong into the household dustbin.  They can be deposited at special refuse centres for used batteries.

Contact your townhall for information.

Taking your Blood Pressure

- Press the -button.
- All digital elements and function symbols light up momentarily to indicate that the unit is functioning correctly.

Placing the Cuff:

- Remove any tight clothing before taking the blood pressure to help ensure that there is no pressure on the arm that may restrict the normal blood circulation.

- Open cuff into a ring shape and pull onto your upper arm until the lower cuff edge lies 2-3 cm (1 inch) above the elbow. The cuff is positioned correctly if the special marking is located directly on the artery. The metal ring must never come to lie on the artery as this may cause a false reading.

Most persons have a higher blood pressure on their left arm. Therefore the measurement should be taken on the left arm. If the blood pressure is higher on the right arm, the measurement should be taken on the right arm.

- Lightly pull the cuff end which feeds through the metal ring and wrap it over the outside of the upper arm. The velcro fastening is then secured with a light pressure.

Do not draw the cuff too tight. Allow a space of ca 2-finger width between arm and cuff.

- Note on cuff sizes:
This unit must only be used with the following cuff types:

Type	Size of Rubber Bag	Arm Circumference	Order No.
CA01	13 x 25 cm	22 – 32 cm	143-4-762*
CA02	14 x 30 cm	32 – 48 cm	143-4-763

* comes as standard


- The "O" in the display and three bleeps indicate that your BP is now ready for use.
- Place the arm with the cuff on the table, holding it relaxed and at a slight angle.
- Now take the rubber ball into your right hand and inflate the cuff to approx. 180 mmHg. A triple signal tone occurs. The rising cuff pressure is shown on the digital display. In case the pressure in the cuff is not high enough, the symbol ▲ appears during the measurement. Please inflate the cuff by approx. 40 mmHg more. This procedure may be repeated until the pressure in the cuff is higher than the systolic blood pressure (symbol ▲ disappears and the triple signal tone occurs).
- At the following measurement the necessary inflation pressure is adjusted by the **memomatic** function, i.e. inflation pressure = **SYSTOLIC** of last measurement taken +40 mmHg). Please inflate until the triple signal tone indicated the correct pressure.

- **Important:** The additional increase in the cuff pressure can cause faulty readings. If this occurs, repeat the measurement after a 2 minute interval and this time inflate straight to above your upper (systolic) blood pressure value.
- After inflation the air is automatically released. The decreasing cuff pressure is registered on the left side of the display.
- As soon as the first pulsations are registered, the symbol  blinks in rhythm with the pulse.

Keep perfectly still and do not speak while taking your blood pressure.

- As soon as the registered pulse is ready to be evaluated by the microprocessor the unit emits a bleep, again in pulse rhythm. At this point it is crucial that the arm with the cuff is not being moved.
- The completion of the measurement is indicated by an extended bleep while the two measurement values in the display remain static.
- Symbol  lights up to indicate that the quick air release can now be initiated. To do this, press the valve button until all the air has been expelled from the cuff.
- If symbol "" appears following a completed measurement, it is recommended that the measurement is repeated while holding the arm motionless. A renewed display of symbol "" following the repeat measurement with the arm being held motionless, can be an indication of irregular heart beat. Discuss this matter with your GP on your next visit.

In the display appear your BP values:

- **SYSTOLE** (upper value)
– left display
- **DIASTOLE** (lower value)
– right display

- **PULSE** (per minute)
– automatically alternating with the above display

- If no further measurement is required, remove the cuff from your arm.
- Enter the obtained values in the accompanying record card. The values will stay in the display for 1 minute after which time the unit switches itself off.
- To repeat a measurement, press the -button again.
- After an interval of **two minutes** a repeat measurement can take place. Only the cuff needs to be inflated again.

Interruption of Measurement

To interrupt a measurement you can press the button of the release valve until all air has been expelled from the cuff. (To repeat a measurement, simply reinflate the cuff.)

Error Messages

If the unit does not emit any rhythmic pulse bleeps it means that the registered pulse is too weak and there will be no measurement results.

The symbols *Err* or *Err CUF* in the display indicate that no blood pressure values or pulse were available for measurement.

- Tighten the cuff around the arm and take a new measurement.

“♥” : Irregular puls or incorrect taking of blood pressure (i.e. movement during measurement)
→ Repeat measurement holding the arm absolutely still. If the “♥” symbol is displayed again, it can be that an irregular heart beat is present. Consult your GP.

Cleaning of Unit and Cuff

- To clean the unit only use a soft dry cloth.
- Small stains on the cuff can be removed with any proprietary cleaner.

Guarantee information, Care of the Environment

For details on the conditions of guarantee, please refer to the guarantee card supplied.

A claim for guarantee can only be allowed if the Certificate of Warranty, completed and stamped by the dealer, is enclosed with the instrument.

For both guarantee and repair work, please send the instrument carefully packed and postage paid either to your authorized dealer or directly to:

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Used batteries and defective electronic devices do not belong into the household waste.

They can be deposited at special centres. Please contact your local authorities for information.

Technical Data

Measuring Principle:

Oscillometric

Measurement range:

40 to 280 mmHg; 40 to 200 Pulse/min.

Accuracy of

pressure indicated ± 3 mmHg, pulse $\pm 5\%$

Display:

LCD

Working Conditions:

Air temperature: 10 – 40 degrees Celsius

Humidity: 10 – 85%

Power:

DC 1,5 V

battery 1x1,5 V

IEC R 6 or LR 6

Battery Control:

Symbol in Display

Weight:

0,1 kg

Classification:

Type BF ()

Clinical Test:

In accordance with DIN 58130,

Test procedure N 2: "Comparative BP measurements at left and right arm"

All rights to technical modifications reserved.

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment.

The use of accessories and cables other than those specified (other than boso original parts) may result in increased emissions or decreased immunity of the unit.

Guidance and manufacturer's declaration – electromagnetic emissions

The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The boso unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The boso unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	n.a.	
Voltage fluctuations/flicker emissions IEC 61000-3-3		

Recommended separation distances between portable and mobile RF communications equipment and the boso unit

The boso unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the boso unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the boso unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
W	$d = 1,2 \sqrt{P}$	$d = 1,2 \sqrt{P}$	$d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitter rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.

Immunity test	IEC 60601-test level	Compliance level	Electromagnetic environment – guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any party of the boso unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p>
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
			<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the boso unit is used exceeds the applicable RF compliance level above, the boso unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the boso unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity

The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.

Immunity test	IEC 60601-test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	n.a.	
Surge IEC 61000-4-5	± 1 kV differential mode ±2 kV common mode	n.a.	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (> 95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	n.a.	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	

NOTE : U_T is the a.c. mains voltage prior to application of the test level.

A) Function testing

A function test can only be carried out on a person or with a simulator.

B) Testing of Density of Pressure Circuit and Divergence from Pressure Indication

1. Remove batteries.
2. Hold -button pressed and at the same time insert batteries.
3. Release -button. The unit now is ready for testing. The displays of SYS and DIA now show the actual pressure.
4. Prior to testing the density, open the battery compartment and take the T-connector off the unit. Wrap the cuff around a fixed cylinder and connect it with the T-connector (without deflation valve) to the pump ball and the unit. Mind the setting time for cuff: at least 30 seconds. Accomplish testing in the usual manner.

5. Prior to testing the divergence of the pressure indication, open the battery compartment and take the T-connector off the unit. Connect the unit with the normal unit and a pump ball. Accomplish testing in the usual manner.

6. Switch off unit.

7. Re-connect the cuff and the pumpball with the T-connector.

C) Placing of Safety Stamps

For complete safety cover up the drill hole for the base screw.

Alternatively, connect upper housing and lower housing with a safety stamp placed directly on one of the housing edges.

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BOSCH + SOHN GmbH u. Co. KG