(D) Gebrauchsanweisung

(GB) User Instructions

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(E) Manual del Usuario

boso
carat
professional

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Package

- 1 Blood Pressure Unit
boso carat professional
- 3 Cuffs CA01, CA02 and CA03
- 1 Guarantee Certificate
- 4 Batteries LR 6
- 1 Power Supply Unit
- 1 User Instructions

Preliminary Notes

Your **boso carat professional** unit is a state of the art blood pressure computer. It works on the **oscillometric measuring principle**. The oscillations created by the pulse and transmitted via the cuff are stored and evaluated by a micro processor. Therefore there is no need for a microphone, the positioning of which is vital if one is to get reliable readings.

The boso carat professional conforms to the requirements of the European Directives based on the Law governing Medical Products (CE mark), and to the European Standard EN 1060, Part 1, on „Non-Invasive Blood Pressure Instruments – General Requirements“ and Part 3 „Additional Requirements for Electronic Blood Pressure Systems.“

Regular checks on the accuracy of the measurement system should be carried out in accordance with the valid directives of the individual countries, governing this issue.

The instructions for re-calibration checks are on page 32 of this booklet.

Important Observations to be followed when Self Testing

- The measurement is taken on the free upper arm. Please ensure that there is no restriction from a tight, pushed-up sleeve. It may be necessary to remove the garment.
- During measurement the arm has to be kept absolutely still.
- Between successive measurements an interval of at least two minutes must be allowed.
- Irregular heart beats can influence the sensitivity of the unit and result in inaccurate measurement values being displayed. Adverse results can also be experienced by patients with a pacemaker who have a weak pulse. The BP unit itself has however no detrimental effect on the pacemaker.
- Using your blood pressure monitor in the vicinity of strong electro-magnetic fields, (ie mobile telephones) can also cause faulty readings (page 33-35).

Starting up

Fitting the Batteries

The battery compartment is located on the underside of the unit. Insert batteries as shown in the illustration.

Please note: If batteries are inserted incorrectly the unit does not function and this can lead to warming and leakage.

- Only use high quality leak proof batteries of the correct specifications (see Section Technical Data) and never mix used with new batteries or batteries of different manufacture.
- Remove batteries if the unit is not to be used for any length of time.
- The unit shows the following function symbols:
 - ☐  battery fully loaded
 - ☐  battery partially loaded
 - ☐  replace battery soon
 - ☐  flashing! No further measurements possible; replace battery.

● Care for your environment

Used batteries and recharging units do not belong into the household dustbin.

They can be deposited at special refuse centres for used batteries. Contact your townhall for information.

● Change batteries only when the unit is switched off!

The measurement values will be deleted if the batteries are not replaced within 30 seconds. With the removal of the batteries while the unit is switched on, the memory store is erased immediately.

How to take your own Blood Pressure

Fitting the cuff

- Open cuff into a ring shape and pull onto your upper arm until the lower cuff edge lies 2-3 cm (1 inch) above the elbow. The cuff is positioned correctly if the special marking is located directly on the artery. The metal ring must never come to lie on the artery as this may cause a false reading.

Most persons have a higher blood pressure on their left arm. Therefore the measurement should be taken on the left arm. If the blood pressure is higher on the right arm, the measurement should be taken on the right arm.


- Gently pull the cuff end that feeds through the metal ring and place it up and outside around your arm. Close velcro fastening with a light pressure. The cuff must not be too tight. There should roughly be space for two fingers between arm and cuff.
- Note on cuff sizes:
This unit must only be used with the following cuff types:

Type	Size of Arm Rubber Bag	Circumference	Order Number
CA01	13 x 25 cm	22 – 32 cm	143-4-766*
CA02	14 x 30 cm	32 – 48 cm	143-4-767*
CA03	9 x 18 cm	16 – 22 cm	141-4-772*

* comes as standard


Important :

A measurement should not take longer than 2 minutes. In case of an error the measurement can be interrupted at any time by opening the velcro fastening of the cuff.


Taking your Blood Pressure

- Attach cuff the BP monitor by fitting the plug firmly into (click) the opening at the left side of the unit.
- Place the arm with the cuff on the table holding it relaxed and slightly bent so that the cuff is line with your heart.

- Press the START-button. All digits and function symbols in the display are momentarily lit up to indicate that the unit is now ready for use.

- The unit has an „intelligent“ automatic system for inflating the cuff to the correct pressure. The rising cuff pressure is displayed by the changing digits in the window.

- When the correct pressure is reached the pump switches off and the air is automatically released from the cuff.


- While the blood pressure is being taken the drop in cuff pressure is again digitally displayed and symbol ♥ is lit up. Remember to remain perfectly still and not to speak.

Display of Measured Values

- As soon as the first pulsations are registered, symbol  flashes in Pulse rhythm.
- After completion of a measurement the built-in magnetic valve opens automatically to allow rapid deflation of the cuff.
- The latest measurement is automatically retained in the memory.
- Invalid measurements are not stored in the memory.
- If symbol  appears following a completed measurement, it is recommended that the measurement is repeated while holding the arm motionless.
A renewed display of symbol  following the repeat measurement with the arm being held motionless, can be an indication of irregular heart beat.
- The unit switches off automatically after approximately one minute. To repeat the procedure, once more press START.
- If no further measurement is required, remove cuff.

Interruption of a Measurement :

A measurement can be interrupted at any stage by simply pressing the START-button. This triggers the automatic air release from the cuff. For further measurements press START again.

Recalling the last recorded measurement

- To recall the result of the last measurement, press the START-button for approx. 3 seconds, while the unit is switched off.
- The result of the last measurement appears for approx. 3 seconds. Afterwards the unit switches off automatically.
- To start a new measurement, simply press „START“.
- **Change batteries only when the unit is switched off!**
The measurements values will be deleted if the batteries are not replaced within 30 seconds. With the removal of the batteries while the unit is switched on, the memory store is erased immediately.


Mains Operation

- A mains socket connector is located on the rear side of the unit. For operating the unit from the mains use only the **boso power supply pack (order number 410-7-153)**. This power supply pack produces a rectified output of the correct polarity. Other commercially available power supply units can cause damage to the electronic components resulting in the guarantee becoming void.

Important: If no batteries are in the unit, the memory store is cleared as soon as the power supply lead is disconnected. Insert batteries in order to save the memory store.
Switch off the unit, disconnect the connector from the unit and then the adapter from the wall socket.

Error Messages

- Should there be any problem while taking your blood pressure that impedes correct monitoring, the display will show an error message instead of a BP reading.

- Meaning of the error messages

Err: cuff pressure varied while blood pressure is taken
 → hold arm still

– or –

no valid pulse readings could be obtained
→ check position of cuff

– or –

Systolic – Diastolic pressure ≤ 10 mm/g
→ check position of cuff

Err Incorrect inflation

CUF: → cuff possibly not firmly enough positioned

“♡” : Irregular pulse or incorrect taking of blood pressure (i.e. movement during measurement)
→ Repeat measurement holding the arm absolutely still. If the “♡” symbol is displayed again, it can be that an irregular heart beat is present.

Cleaning of BP Unit and Cuff

- To clean your unit use a soft dry cloth only.
- Small stains can be removed with a proprietary stain remover.

Guarantee information, Care of the Environment

For details on the conditions of guarantee, please refer to the guarantee card supplied.

A claim for guarantee can only be allowed if the Certificate of Warranty, completed and stamped by the dealer, is enclosed with the instrument.

For both guarantee and repair work, please send the instrument carefully packed and postage paid either to your authorized dealer or directly to:

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Used batteries and defective electronic devices do not belong into the household waste.

They can be deposited at special centres. Please contact your local authorities for information.

Technical Data:

Measurement system: Oscillometric

Measurement range: 40 – 240 mmHg, 40 – 180 Puls/min.

Cuff pressure: 0 – 320 mmHg

Memory store: 1 measurement

Display: LCD

Working conditions:

Environmental temperature: Room temperature 10 – 40°C

Relative humidity 10 – 85 %

Storing conditions:

Environmental temperature -5°C to +50°C

Relative humidity 85 % maximum

Power supply:

DC 6 V (batteries 4 x 1,5 Mignon IEC LR 6, alkaline manganese)

Alternative special option: Main unit DC 6 V, Order-No. 410-7-153 (Polarity: outside MINUS, inside PLUS) 

Battery check: Symbol display in LCD window

Weight: 600 g without batteries

Dimensions (B x H x T):

150 mm x 160 mm x 120 mm

Classification: Class II ()

Typ BF ()

Clinical Test (DIN 58130):

Accuracy complies with EN 1060 part 3

Maximum deviation of cuff pressure: ± 3 mmHg

Maximum deviation of pulse rate: ± 5 %

The right to alter technical specifications is reserved.

A) Function testing

A function test can only be carried out on a person or with a suitable simulator.

B) Testing of the Density of Pressure Circuit and Divergence in the Pressure Indication

Please note:

- a) When in the measurement mode pressure is increased to above 320 mmHg, the quick release valve is triggered and the pressure circuit opened. When during testing pressure is increased to above 320 mmHg the measurement display will flash.
- b) Prior to testing the connector piece has to be removed from the cuff tube. Then the short connecting head of the connector plug – up to now attached to the unit and its corresponding socket at the unit – has to be attached to the cuff tube, and the long connecting head of the connector piece to be inserted into the unit.
Lastly the pump ball is to be fitted into the pressure circuit.

Testing

- 1.) Remove batteries.
- 2.) Hold START-button pressed down and insert batteries.
- 3.) Release START-button.
The unit is now ready for testing. The displays of SYS, DIA and PUL now show the actual pressure.
- 4.) Testing for divergence of Pressure indication and Density of Pressure circuit (setting time for cuff at least 30 seconds) can now be carried out in the usual manner.
- 5.) Switch off unit and re-establish all normal connections.

C) Safety Seal

As a measure to safeguard the correct calibration of the BP unit the upper and lower parts of the housing can be connected with a safety seal. Alternatively, a hole drilled for the connecting screw at the back of the unit can be sealed for safety.

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according the EMC information provided in the following.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment.

The use of accessories and cables other than those specified (other than bosco original parts) may result in increased emissions or decreased immunity of the unit.

Guidance and manufacturer's declaration – electromagnetic emissions		
The bosco unit is intended for use in the electromagnetic environment specified below. The customer of the user of the bosco unit should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The bosco unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The bosco unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the bosco unit			
The bosco unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the bosco unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the bosco unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitter rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration – electromagnetic immunity

The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the boso unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p>
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{eff}	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
			<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the boso unit is used exceeds the applicable RF compliance level above, the boso unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the boso unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity			
The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	n.a.	
Surge IEC 61000-4-5	± 1 kV differential mode ±2 kV common mode	± 1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (> 95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	< 5% U_T (> 95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the boso unit requires continued operation during power mains interruptions, it is recommended that the boso unit has to be powered from a uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	
NOTE : U_T is the a.c. mains voltage prior to application of the test level.			

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BOSCH + SOHN GmbH u. Co. KG