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ⒺN User Instructions

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**Blutdruckmessgerät mit Stethoskop**  
**Sphygmomanometer with Stethoscope**  
**Tensiomètre avec stéthoscope**  
**Sfigmomanometro con stetoscopio**  
**Tensiómetro con fonendoscopio**

**boso** med I  
BS 90  
varius privat

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## Contents of the kit

- 1 blood pressure instrument with cuff and stethoscope
- 1 set of instructions for use
- 1 case

## Symbols on the blood pressure monitor

Symbol	Function/meaning
	Device complies with the European Medical Devices Directive
	Year of manufacture

## Preliminary remarks

This blood pressure instrument complies with the international standard ISO 81060-1.

All Calibration testing – at least every two years –, has to be carried out either by the manufacturers or the local authorised dealer in accordance with the relevant distributor directives.

A modification of the device and/or the accessories is not permitted. This can lead to measurement errors.

The device can be used by any user skilled in auscultatory blood pressure measurement.

### **Purpose:**

Non-invasive recording of systolic and diastolic blood pressure in humans. The instrument can be used with any arm circumference as shown on the corresponding cuff.

### **Important**

Keep this unit from being knocked or dropped, and protect it from dirt and moisture. Care has to be taken not to puncture the cuff with any sharp instruments (Scissors, Needles).

---

Do not inflate to above 300 mmHg !

Measurement time should be ca. 2 minutes maximum. Between measurements allow an interval of at least 2 minutes.

**Please use only boso cuffs.**

The cuff must be selected such that it fits the arm circumference.

**What you should know about your blood pressure**

Your heart is your engine for life. You know that. But do you also know how it works ?

Constant overcharge means increased wear and tear, or in other words, a shorter life expectancy.

Your blood pressure tells you the power your heart produces, or has to produce.

A constantly raised blood pressure means more stress, a lower blood pressure, less stress on your heart.

Blood pressure which is constantly above the normal rates leads – in time – to serious health problems, risks such as arteriosclerosis, heart attacks, strokes, and sometimes kidney failure, and retinitis.

Early diagnosis and treatment of lasting hypertension

(high blood pressure) can help prevent these risks.

It is advisable therefore – especially from the age of 40 – to have your blood pressure checked several times during the year. If high blood pressure is recognised, you should monitor your blood pressure for several consecutive days, mornings and evenings, at roughly the same time and in the same body position (sitting or lying down), generally after ca. 5 minutes rest. This ensures your circulation is normal, this providing equal conditions for comparable measurements.

If your blood pressure remains constantly above normal over a period of several days, you should not hesitate to consult your doctor and show him your record card with the measurements taken.

He will make the correct diagnosis and if necessary, prescribe the necessary medication to lower your blood pressure.

Do not get alarmed if on occasions you measure a higher blood pressure than normal. Blood pressure changes constantly, because the heart adapts to the physical demands of the body.

It is perfectly normal for blood pressure to increase substantially after physical exercise, as well as after exceptional psychological influences such as shock, fear, excitement etc.

Only if a constantly high blood pressure is found in complete rest, you have real hypertension, which further increases dangerously under physical or mental stress.

Although hypotension (low blood pressure) can cause uncomfortable symptoms like vertigo and general exhaustion, the sub-normal blood pressure does not represent a serious health risk. It can generally be regarded as a „resting mode” for your engine heart, with a higher life expectancy.

The World Health Organization (WHO) has set the following guideline for the assessing of blood pressure values:

	Systolic	Diastolic
Too high	over 140 mmHg	over 90 mmHg
Normal – borderline	130 to 139 mmHg	85 to 89 mmHg
Normal	120 to 129 mmHg	80 to 84 mmHg
Optimum	up to 119 mmHg	up to 79 mmHg

Whether medical treatment is necessary depends not only on the blood pressure, but also on the patient’s risk profile. Please contact your GP if one of the values (systole, diastole) consistently exceeds the cut-off point for high blood pressure.


Self-measurement of blood pressure often produces somewhat lower values at home than at the doctor’s. The German Hypertension League has therefore defined lower cut-off values for home measurements:

Home measurement: 135/85 mmHg

Measurement in the doctor’s surgery: 140/90 mmHg

## How to measure your Blood Pressure

- 1.) Sit down comfortably and relaxed in a chair with arm rests, or at a table onto which you can place your arm at a slight angle.
- 2.) not cross his/her legs
- 3.) rest his/her back and arms on the chair or another surface
- 4.) have the centre of the cuff placed on the upper arm at the level of the heart
- 5.) Uncover your left arm and place your forearm on the table. On the inside of the upper arm, ca. 2,5 cm above your elbow, your can feel your pulse (aorta) with your finger tips.
- 6.) Make sure that there is no pressure on the arm from the pushed-up sleeve to restrict the circulation. Tight clothing is best removed before taking your blood pressure.

- 7.) Slide the open cuff over your upper arm until the lower cuff edge is ca. 2 – 3 cm away from the elbow. The rubber tubes should extend down from the cuff. The special care that the membrane in cuff is positioned directly over the pulse located earlier. The more precise you position the membrane the clearer and louder you will be able to hear the sounds of the artery.
- 8.) relax as much as possible and keep silent while his/her blood pressure is being measured

9.) be calm for 5 minutes before the first reading

The following phases are recommended in the case of auscultatory measurement


- adults: phase V (K5)
- children aged 3 to 12: phase IV (K4)
- pregnant women: phase V (K5), except when sounds are audible when the cuff is deflated: in this case phase IV (K4) Korotkoff sounds should be used.

K5 is the point at which tones heard through the stethoscope are no longer audible

K4 is the point at which tones heard through the stethoscope change from a clear sound to a muffled sound.

**Inscriptions:**

“mmHg” is the unit in which blood pressure is measured  
“shock protected” = the instrument has been tested and found to work correctly after being dropped from a height of up to one metre.

10.) **D-Ring Velcro Cuff**

The cuff must not be fitted over wounds as this could lead to further injury.

Blood pressure must always be measured on the arm with the higher pressure values.

To do this, first measure blood pressure on both arms and subsequently always on the arm which produced the higher reading.

---

Draw the cuff close. Place the free end outwards over the upper arm (see Illustration) and close the velcro fastening.

The cuff should fit firmly around the upper arm, but not exert any pressure on the artery before inflation. If you can insert a finger between cuff and upper arm, the cuff sits correctly.

- 11.) Put on the stethoscope allowing the earpieces of the head frame to fit well inside your ears. These arms can be rotated within their spring so that they can be positioned individually.
- 12.) Now place the arm with the cuff back on the table, comfortably, slightly angled, resting. Take the manometer into your right hand, holding it around the pump, with the thumb and forefinger on the valve screw. Close the valve with a right turn and pump up the cuff steadily, watching the manometer.  
Continue to pump until the indicator on the manometer exceeds your usual blood pressure by ca. 30 mmHg.

Should you still hear pulsing sounds through the stethoscope, continue to inflate until these sounds stop, which indicates that the circulation is disrupted.

- 13.) Open the valve very slightly with thumb and forefinger and watch the slow decrease in pressure on the manometer. The rate of deflation between systole and diastole should be in the range of 2 – 3 mmHg per second, ie there should be an interval of ca. 20 – 25 seconds for the pointer on the scale to move from 150 to 100. To avoid a reading error, read the pressure values perpendicular to the scale.

- **When you hear the first clear arterial sound (knocking sound) through the stethoscope you can read the systolic (upper) blood pressure on the manometer scale.**
- **When the arterial sound stops, you can read the diastolic (lower) blood pressure.**

- 14.) After you have taken your readings, open the valve screw fully turning it to the left, to allow all air to be released swiftly from the cuff.  
If you want to repeat your measurement, allow at least two minutes for your circulation to get back to normal. If no further measurements are required, remove the cuff.

## Technical Data

Measurement range: 0 – 300 mmHg

Accuracy of Pressure shown:  $\pm 3$  mmHg or 2% of the reading (whichever is greater)

Storing conditions:

$-30^{\circ}\text{C}$  –  $70^{\circ}\text{C}$ , 15 – 85% rel. humidity

Operating conditions:

$10^{\circ}\text{C}$  –  $40^{\circ}\text{C}$ , 15 – 85% rel. humidity

## Cleaning of Unit and Cuff

To clean the unit only use a dry, soft cloth. To clean the cuff – remove the rubber bag. The exterior of the cuff can then be washed by hand at a maximum temperature of 30 degrees Celsius.

## Disinfection

For disinfectant wipes (at least 5 minutes exposure time) of the device and the cuff, we recommend the disinfectant mikrozyd sensitive liquid (Schülke & Mayr).  
To disinfect the cuff, we recommend spray disinfection.

---

## Guarantee information

We give 3 years warranty from the date of purchase. The purchase date has to be proven by the invoice. Within the warranty period defects are eliminated free of charge. After repairs the warranty period is not extended on the whole unit, but only to the replaced components.

Excluded from the warranty are parts subject to normal wear and tear (e.g. cuff), transport damages and any damage caused by improper handling (e.g. non-compliance with the instructions for use). Damages due to disassembly by unauthorized persons are also excluded from warranty.

No claims for damages against us are substantiated by the warranty.

In the case of justified warranty claims the device has to be sent along with the original invoice to:

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## Disposal

Please ensure that the device is disposed of according to all regional and national environmental regulations.

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**BOSCH + SOHN GmbH u. Co. KG**