
boso mercurius E

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(DE)

2

User
Instructions

(EN)

36

Mode
d'emploi

(FR)

70

Käyttöohje

(FI)

104

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Table of contents

EN


Product features	38
Contents of package	39
Symbols of the blood pressure monitor	40
Quick guide	44
Purpose / Contraindications	45
Introduction	46
Starting up	48
Inserting batteries	48
Battery status	49
General instructions for measurement	50
Preparing to measure your blood pressure	51
Note on cuff sizes	51

Table of contents

37

Measuring your blood pressure.....	52
Attaching the cuff.....	53
Pulse display.....	58
Error messages	60
Cleaning and Disinfection.....	61
EN	
Warranty Conditions	62
Customer service.....	63
Disposal.....	64
Technical data.....	65
EMC-Notes.....	67

(EN)

Contents of package

39


- 1 **boso mercurius E**
blood pressure monitor
+
1 standard cuff
with pump ball

2 batteries LR 6 (AA)

1 User
instructions

EN

40 Symbols on the blood pressure monitor

Symbol	Function/meaning
	turn device on/off
	Manual marking of systole and diastole in the display
	Serial number
	Year of manufacture
	Type BF: The device, cuff and tubing have been designed to protect the patient against electric shock.

Symbols on the blood pressure monitor

41

Symbol	Function/meaning
	Battery fully charged.
	Battery partly discharged. (EN)
	Replace battery as soon as possible.
	Flashing! No further reading possible. Replace battery.
	Protect against liquids

42 Symbols on the blood pressure monitor

Symbol	Function/meaning
	Pulse not recognized correctly
	Rest of air in system
	Measurement time > 3 min
	Unit defective
mmHg	Cuff pressure unit
bpm	Pulse per minute

Symbols on the blood pressure monitor

43

Symbol	Function/meaning
CE 0124	Device complies with the European Medical Devices Directive.
 	Device must not be discarded with household waste. 
	Manufacturer
	Read instructions for use

Quick guide

- Insert batteries, ensure correct positioning. (P. 48)
- Determine by auscultation the systolic and diastolic values. (S. 57)
- Press  button. (P. 52)
- Attach cuff. (P. 53)
- After the measurement is completed, fully open the air valve (by turning the air valve anti-clockwise) to achieve fast deflation of the cuff. (P. 58)
- Close the air valve next to the pumping ball by turning it clockwise. (P. 57)
- The pulse is shown in the display.
- Inflate the cuff to above the systolic pressure. The air release speed can be adjusted by rotating the valve screw. (P. 57)

EN

⚠ This quick guide does not replace the detailed information on the use and safety of the blood pressure monitor featured in the user instructions.

Please be sure therefore to read the user instructions!

Purpose

Non-invasive recording of systolic and diastolic blood pressure and pulse frequency in humans.

The instrument can be used with any arm circumference as shown on the corresponding cuff.

Contraindications

(EN)

Do not use the device on patients

- who are connected to a heart-lung machine
- who are in a critical medical condition or on an intensive-care ward.

EN

Dear Customer,

We are delighted that you have decided to purchase a boso blood pressure monitor. The boso brand is synonymous with optimum quality and precision and is also number 1 among professionals: 96% of all German general practitioners, physicians and internists work in practice with blood pressure instruments from boso (API survey conducted by GfK 01/2016).

This instrument has also passed our strict quality control procedures and is your reliable partner for monitoring your blood pressure values.

Please read these user instructions before using the instrument for the first time, as correct blood pressure readings can only be obtained if the instrument is operated correctly. In these user instructions, „→“ means an action by the user.

To get help concerning usage or maintenance please contact your dealer or the manufacturer (contact details can be found on the back cover of this manual).

The use of accessories not referred to in these instructions for use can undermine safety.

An unexpected operating condition or an incidence which has worsened the health condition or could have worsened the health condition should be reported to the manufacturer immediately.

If the instrument is to be sold, ensure that these user instructions are enclosed.

If the monitor is used in medical practice (in accordance with the legislation governing the operation of medical devices), it should be calibrated at regular intervals.

(EN)

The boso mercurius E blood pressure monitor is suitable for patients of all ages but not for neonatals.

EN


Illustration 1

- Use only high quality leak-proof batteries of the correct specifications (see "Technical data" P. 65).

⚠ If the batteries are not inserted correctly, the monitor will not work and this may cause the batteries to heat up and leak, damaging the monitor.

⚠ Batteries that are short-circuited can become hot and cause burning.

→ **Inserting the batteries**

The battery compartment is located on the inner side of the monitor. Insert batteries in the battery compartment as shown in illustration 1.

- Remove the batteries if the monitor is not to be used for any length of time.

The monitor has a battery status display

Batteries fully charged.

Batteries partially charged.

Replace batteries soon.

Flashing! No further measurements possible.
Replace batteries.

Care for the environment

Do not dispose of used batteries and charger units in domestic waste.

You can take these to collection sites for used batteries or special waste.

Contact your local authority for information.

EN

50 General instructions for measurement

Always measure your blood pressure from the arm with the higher blood pressure values (where not required otherwise (see page 53)).

(EN)

To do this, measure the blood pressure first in both arms and thereafter always from the arm with the higher blood pressure.

Using the monitor in the vicinity of strong electromagnetic fields (e.g. radiation devices, mobile telephones) can cause it to malfunction. Medical electrical apparatuses are subject to specific precautions as regards electromagnetic compatibility (EMC) (see page 67-69).

Note on cuff sizes:

⚠ The monitor may only be used with original boso cuffs. These should be chosen according to the upper arm circumference printed on them.

Observe the following instructions when measuring your blood pressure:

Before the measurement the patient should:

- avoid smoking or drinking coffee for one hour
- relax for 5 minutes

EN

During the measurement the patient should:

- sit comfortably
- not cross the legs
- place feet flat on the floor
- support back and arms
- not move

52 Measuring your blood pressure

- Open the instrument by pushing the sliding switch to the right. Fold up the instrument.
- Press the  button .

All function symbols light up momentarily to indicate that the unit is functioning correctly.

Illustration 2

The "0" in the lower display indicates that your instrument is now ready for use.

Illustration 3

Measuring your blood pressure


53

Attaching the cuff

The measurement should be taken from the bare upper arm.

⚠ In the case of tightly fitting outer clothing, ensure that the blood supply to the arm is not cut off when the sleeve is rolled up (if necessary remove the article of clothing before taking the measurement).

→ Place the cuff so that the bottom edge is about 2-3 cm above the elbow. The cuff must be positioned so that the marking is over the artery (see illustration 4).

EN

Illustration 4

⚠ Do not attach the cuff over wounds as this can cause further injuries.

⚠ Ensure that the cuff is not attached to an arm where arteries or veins are undergoing or have undergone medical treatment (e.g. shunt).

54 Measuring your blood pressure

(EN)

⚠ For women who have had a breast amputated, do not attach the cuff to the arm on the amputated side of the body.

⚠ During the measurement, medical apparatus that are being used on the same arm at the same time, may malfunction.

⚠ The cuff must not be too tight. There should be enough space to fit roughly two fingers between arm and cuff.

► The cuff should be placed at heart level.


Illustration 5

Do not talk during the measurement.

⚠ Check that the air tube is not kinked during the measurement. This could cause congestion of the blood flow and subsequent injury.

Measuring your blood pressure

55

⚠ The blood pressure measurement should not prevent the flow of blood for an unnecessary length of time (> 2 minutes). If the monitor fails to function correctly, remove the cuff from the arm.

⚠ Excessive frequent measurements can cause injury by restricting blood flow.

⚠ The device may not be used by unsupervised children.

⚠ There are small parts that may cause a choking hazard if swallowed by mistake by infants.

⚠ The device must not be used if the patient experiences excessive pain during measurement.

⚠ The device must not be used in an environment with easily flammable gases (e.g. gases used in anaesthesia) This could cause explosions.

EN

⚠ The device must not be used in an environment with high oxygen levels (e.g. an oxygen tent). This could cause a fire or explosion.

56 Measuring your blood pressure

 If the cuff is soiled by blood or other bodily fluids, dispose of it in accordance with local waste disposal regulations in order to prevent cross-infection.

(EN)

⚠ If the ambient temperature changes significantly, leave the device about an hour to adapt to the new ambient temperature before using it.

⚠ The performance of the device can be affected by excessive temperature, humidity or altitude.

→ Close the air valve next to the pumping ball by turning it clockwise.

→ Inflate the cuff to above the systolic pressure. The air release speed can be adjusted by rotating the valve screw. The World Health

Organisation (WHO) recommends a deflation speed of 2 – 3 mmHg per second..

→ As soon as you realize the systolic auscultation press the "MARK"-button to mark the relevant value on the display.

→ In the moment you realize the diastolic pressure, press the "MARK"-button again to mark the relevant value on the column (maximum 4 marks for diastolic pressure possible)

(EN)

Pulse display

- After the measurement is completed, fully open the air valve (by turning the air valve anti-clockwise) to achieve fast deflation of the cuff.

EN

The pulse is shown in the lower part of the display as well as the marked pressure values that have been marked during measurement.

Illustration 6

If no further measurement is required, switch off the instrument by pressing the button.

To store the instrument, put the cuff back into the compartment and close it.

 Make sure that your fingers are not crushed when closing the device.

(EN)

60 Error messages

If any problems occur during the measurement, the display will show an error message instead of a blood pressure reading.

Meaning of error messages:

(EN)

Cause of Err message:	Troubleshooting:
	pulse not recognized correctly repeat measurement
	rest of air in system please deflate completely
	measurement time > 3 min repeat measurement
	unit defective please sent unit to the service adress

Cleaning Monitor:

Use a soft, dry cloth to clean your monitor.

Cleaning Cuff:

Small stains on the cuff can be removed carefully with a proprietary washing-up liquid.

Disinfection Cuff:

For disinfectant wipes (at least 5 minutes exposure time) of the device and the cuff, we recommend the disinfectant antifect liquid (Schülke & Mayr). To disinfect the cuff, we recommend spray disinfection.

In particular, if the device is used on multiple patients is to pay attention to regular cleaning and disinfections of the cuff.

(EN)

62 Warranty Conditions

(EN)

We give 2 years warranty from the date of purchase. The purchase date has to be proven by the invoice. Within the warranty period defects are eliminated free of charge. After repairs the warranty period is not extended on the whole unit, but only to the replaced components.

Excluded from the warranty are parts subject to normal wear and tear (e.g. cuff), transport damages and any damage caused by improper handling (e.g. non-compliance with the instructions for use). Damages due to disassembly by unauthorized persons are also excluded from warranty.

No claims for damages against us are substantiated by the warranty.

In the case of justified warranty claims the device has to be sent along with the original invoice to:

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Warranty and repair work must be undertaken by trained and authorised personnel. Do not modify this equipment without authorization of the manufacturer.

Send the monitor, carefully packaged and with sufficient postage, to your authorised dealer or directly to:

BOSCH + SOHN GmbH u. Co. KG

Serviceabteilung

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72417 Jungingen

GERMANY

T: +49 (0) 74 77 / 92 75-0

The device must not undergo service or maintenance while it is in use.

EN

Disposal

(EN)

 Used internal and external batteries must not be discarded in household waste. You can take them to a collection point for used batteries or dispose of them in specialist waste disposal facilities. Please contact your local council for more information.

This device is covered by the scope of EC Directive 2002/96/ EG (WEEE). It cannot be disposed of via local authority waste disposal facilities for used electrical equipment. boso has approved a company which will dispose of this device in accordance with the law.

See the back cover of this manual for the address to contact for more information.

Do not throw the packaging material away, send it for recycling.

Anticipated life of the device:
5 years (when used 6 times a day)

Anticipated life of cuff and pump ball:
2 years (when used 6 times a day)

Measurement principle: auscultatory

Measurement range: 20 to 280 mmHg, 40 to 180 pulses/min.

Display: LCD

Operating conditions: environmental temperature +10°C to +40°C
Relative humidity 15 to 85%
Air pressure 800 hPa to 1060 hPa

**Transport/
Storage conditions:** environmental temperature –20°C to +60°C
Relative humidity 10 to 95%

Power supply: DC 3 V (2 x 1.5 V Mignon IEC LR 6 alkaline manganese batteries)

Typical battery life: 2.000 measurement cycles (depending on frequency of use and ambient temperature).

66 Technical data

Battery check: symbol display in LCD window

Weight: 520 g without batteries

 **Dimensions (WxHxD):** 98 mm x 324 mm x 67 mm

Classification: Type BF (██)

Accuracy of pressure indicated: ± 3 mmHg or 2% of the reading
(whichever is greater)

Accuracy of pulse: ± 5%

Applicable standards: ISO 81060-1

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according the EMC information provided in the following.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment.

The use of accessories and cables other than those specified (other than boso original parts) may result in increased emissions or decreased immunity of the unit.

Guidance and manufacturer's declaration – electromagnetic emissions		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The boso unit uses RF energy only for its internal function. There-fore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The boso unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	n.a.	
Voltage fluctuations/flicker emissions IEC 61000-3-3		

Recommended separation distances between portable and mobile RF communications equipment and the boso unit			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitter rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.

Immunity test	IEC 60601-test level	Compliance level	Electromagnetic environment – guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the boso unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p>
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the boso unit is used exceeds the applicable RF compliance level above, the boso unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the boso unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity

The boso unit is intended for use in the electromagnetic environment specified below. The customer or the user of the boso unit should assure that it is used in such an environment.

Immunity test	IEC 60601-test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	n.a.	
Surge IEC 61000-4-5	± 1 kV differential mode ±2 kV common mode	n.a.	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U _T (> 95% dip in U _T) for 0,5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles < 5% U _T (> 95% dip in U _T) for 5 s	n.a.	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	

NOTE : U_T is the a.c. mains voltage prior to application of the test level.

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