boso

BG

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42 Instructions for use

Purpose

Non-invasive blood pressure meter to measure systolic and diastolic blood pressure levels in humans.

Calibration checks

Calibration checks must be performed in accordance with the national regulations in force.

EN

N.B.:

The device must be installed and operated in accordance with the information given in these instructions for use.

The blood pressure meter is suitable for patients of any age. The device is not suitable for infants. If you are pregnant, ask your doctor before using the device. While having their blood pressure taken, patients must sit calmly, with their legs uncrossed and their feet flat on the floor. The back and arms must be supported. Patients must rest in this position for 5 minutes before the first reading is taken. The reading must be repeated if it is obviously incorrect. The process of taking blood pressure must not stop the circulation for an unnecessarily long time (more than 2 minutes). If the device malfunctions, remove the cuff from the arm. Excessively frequent readings can impair

circulation and consequently lead to injury. A gap of at least 2 minutes must be left between 2 readings.

The cuff must not be placed on top of wounds, as this could lead to further injuries. Ensure that the cuff is not placed on an arm of which the arteries or veins are undergoing or have undergone medical treatment, e.g. a shunt. In the case of women who have had a mastectomy, the cuff should not be placed on the arm on the side from which the breast has been removed. Medical devices being used on the same arm at the same time may malfunction.

Blood pressure is a dynamic parameter and can be influenced by various factors: the patient's posture (seated, standing, lying down), movement before or during the reading, physical state (stress, illness, etc.)

This device must be maintained by trained and approved staff. The device must not be modified without the manufacturer's approval.

Please contact your specialist supplier or the manufacturer for help in preparing for first use, using the device or maintenance. Notify the manufacturer immediately of any unexpected operational

44 Instructions for use

status or any event which has or could have caused damage to health.

Wireless communication devices such as wireless home network devices, mobile phones, cordless phones and their base stations and walkie-talkies can affect this blood pressure meter. They should therefore be kept at least 3.3 metres away from the device.

EN

Explanation of symbols:

Battery status display

Reading in progress

Cardiac dysrhythmia during reading

Stored value

Introduction

Two values need to be measured to determine blood pressure. **The systolic (higher) blood pressure.** This is produced when the heart contracts and pushes blood into the blood vessels. **The diastolic (lower) blood pressure.** This is produced when the heart muscle is relaxed and is refilling with blood.

Blood pressure readings are expressed in mmHg (mm mercury column).

The following thresholds have been set by the World Health Organization (WHO). Ideally, these values should be based on the average of several different readings.

	systolic / mmHg	diastolic / mmHg
clearly elevated	140 or above	90 or above
high end of normal	130-139	85-89
normal	120-129	80-84
optimum	up to 119	up to 79

46 Instructions for use

Some devices have four coloured markings to the left of the display screen, corresponding to the WHO thresholds. The display screen shows both the average value and a symbol next to one of the four coloured markings. For example, a symbol next to the red marking points to clearly elevated blood pressure.

EN

Preparing for a reading

➔ Place the batteries in the battery compartment as shown on the diagram. For upper arm devices: connect the cuff to the device. Please bear in mind that only original boson cuffs of appropriate size may be used (check the cuff imprint).

➔ **General principles:** The reading should be taken with the patient seated and relaxed. In individuals who have higher blood pressure in the left arm, readings should be taken on the left arm or wrist. However, in individuals who have higher blood pressure in the right arm, readings should be taken on the right arm or wrist.

➔ The arm on which the cuff has been placed must not be moved during the reading, and the patient must not speak.

Reading taken on the upper arm:

- ➔ Place the bare arm on the table or armrest and bend it slightly. Open the cuff to form a ring shape and slip it over the upper arm until the lower edge of the cuff reaches a point about 2-3 cm above the elbow. The cuff must be positioned so that the marking lies on the artery of the upper arm.
- ➔ Make sure that the air tube is not kinked during the reading. The resulting blood congestion could cause injury.

- ➔ Press the start button.
- ➔ Wait for the result to be displayed.

Reading taken on the wrist:

- ➔ Position the device on the bare wrist. It must be placed on the inner side of the wrist.
- ➔ Bend the arm. The wrist must be at the level of the heart during the reading.
- ➔ Press the start button.
- ➔ Wait for the result to be displayed.

48 Instructions for use

Arrhythmia detection (optional):

If irregular pulsations occur during the reading, a heart in brackets is displayed along with the result once the reading has been completed. If this symbol appears frequently while the arm is still, this may indicate heart rhythm abnormalities. Please discuss this when you next see your doctor.

EN

Cleaning and disinfecting the device and the cuff:

Use a soft, dry cloth to clean the device. Small spots on the cuff can be removed carefully using ordinary washing-up liquid. We recommend using antiseptic liquid

disinfectant (Schülke & Mayr) to disinfect the device by wiping it, leaving the disinfectant on for at least 5 minutes. We recommend a disinfectant spray to disinfect the velcro on the cuff. Ensure that the cuff is regularly cleaned and disinfected, especially if the device is being used by several people.

Disposal

The device and batteries must not be discarded in household waste. When the device has reached the end of its useful life, it must be sent to an electronic waste collection facility. Please note that batteries and rechargeable batteries must be disposed of separately. (Your local council's waste collection facility).

Anticipated life of the device: 10 years
Anticipated life of cuffs: 10,000 readings
Typical battery life: 500 readings
(depending on inflation pressure and frequency of use)

50 Technical data

Principle of measurement:	Oscillometric
Measurement range:	40 to 240 mmHg, 40 to 200 pulses/min.
Cuff pressure:	0 to 300 mmHg
Display:	LCD
Protection against foreign bodies and water:	IP20
 Operating conditions:	Ambient temperature +10°C to +40°C - relative humidity 15 to 85%
Storage conditions:	Ambient temperature -10°C to +60°C - relative humidity 15 to 85%
Power supply:	DC 3/6 V (batteries 2/4 x 1.5 V AAA/AA IEC LR 03/06, alkaline manganese), optional extra feature: Power pack DC 6 V (order no. 410-7-150)

Cause of Err display	Action:
Cuff pressure fluctuates too much while the reading is being taken	Keep arm still
No measurable pulsations	Check cuff position and repeat reading
Difference between systolic and diastolic pressure too low (sys - dia < 11 mmHg)	Check cuff position and repeat reading
Cause of Err CuF display:	Action:
Error during inflation	Cuff may be too loose, tighten cuff

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BOSCH + SOHN GmbH u. Co. KG