- ⓓ Gebrauchsanweisung
- ⓖⓔ User Instructions
- ⓕ Mode d'emploi

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Introduction

This blood pressure monitor conforms to the European requirements that are the basis of the Medical Products Law (CE), as well as to the European Standard „Europanorm“ EN 1060, Part 1 „Non-Invasive blood pressure monitoring devices – General requirements“ and part 2 „Additional requirements for mechanical blood pressure monitoring devices“.

All Calibration testing – at least every two years –, has to be carried out either by the manufacturers or the local authorised dealer in accordance with the relevant distributor directives.

Important

Keep this unit from being knocked or dropped, and protect it from dirt and moisture. Care has to be taken not to puncture the cuff with any sharp instruments (Scissors, Needles).

Do not inflate to above 300 mmHg !

Measurement time should be ca. 2 minutes maximum. Between measurements allow an interval of at least 2 minutes.

Please use only boso cuffs.

What you should know about your blood pressure

Your heart is your engine for life. You know that. But do you also know how it works ?

Constant overcharge means increased wear and tear, or in other words, a shorter life expectancy.

Your blood pressure tells you the power your heart produces, or has to produce.

A constantly raised blood pressure means more stress, a lower blood pressure, less stress on your heart.

Blood pressure which is constantly above the normal rates leads – in time – to serious health problems, risks such as arteriosclerosis, heart attacks, strokes, and sometimes kidney failure, and retinitis.

Early diagnosis and treatment of lasting hypertension (high blood pressure) can help prevent these risks.

It is advisable therefore – especially from the age of 40 – to have your blood pressure checked several times during the year. If high blood pressure is recognised, you should monitor your blood pressure for several consecutive days, mornings and evenings, at roughly the same time and in the same body position (sitting or

lying down), generally after ca. 3 minutes rest. This ensures your circulation is normal, this providing equal conditions for comparable measurements.

If your blood pressure remains constantly above normal over a period of several days, you should not hesitate to consult your doctor and show him your record card with the measurements taken.

He will make the correct diagnosis and if necessary, prescribe the necessary medication to lower your blood pressure.

Do not get alarmed if on occasions you measure a higher blood pressure than normal. Blood pressure changes constantly, because the heart adapts to the physical demands of the body.

It is perfectly normal for blood pressure to increase substantially after physical exercise, as well as after exceptional psychological influences such as shock, fear, excitement etc.

Only if a constantly high blood pressure is found in complete rest, you have real hypertension, which

further increases dangerously under physical or mental stress.

Although hypotension (low blood pressure) can cause uncomfortable symptoms like vertigo and general exhaustion, the sub-normal blood pressure does not represent a serious health risk. It can generally be regarded as a „resting mode“ for your engine heart, with a higher life expectancy.


The World Health Organization (WHO) has set the following guideline for the assessing of blood pressure values:

	Systolic	Diastolic
Too high	over 140 mmHg	over 90 mmHg
Normal – borderline	130 to 139 mmHg	85 to 89 mmHg
Normal	120 to 129 mmHg	80 to 84 mmHg
Optimum	up to 119 mmHg	up to 79 mmHg

If your doctor diagnoses high blood pressure, you must strictly adhere to the prescribed dosage of medication, and do not change this when you occasionally get lower blood pressure readings.

How to measure your Blood Pressure

- 1.) Sit down comfortably and relaxed in a chair with arm rests, or at a table onto which you can place your arm at a slight angle.
- 2.) Uncover your left arm and place your forearm on the table. On the inside of the upper arm, ca. 2,5 cm above your elbow, you can feel your pulse (aorta) with your finger tips.
- 3.) Make sure that there is no pressure on the arm from the pushed-up sleeve to restrict the circulation. Tight clothing is best removed before taking your blood pressure.
- 4.) Slide the open cuff over your upper arm until the lower cuff edge is ca. 2 – 3 cm away from the el-


bow. The rubber tubes should extend down from the cuff.

The special care that the membrane in cuff is positioned directly over the pulse located earlier. The more precise you position the membrane the clearer and louder you will be able to hear the sounds of the artery.

5.) **D-Ring Velcro Cuff**

Draw the cuff close. Place the free end outwards over the upper arm (see Illustration) and close the velcro fastening.

The cuff should fit firmly around the upper arm, but not exert any pressure on the artery before inflation. If you can insert a finger between cuff and upper arm, the cuff sits correctly.

6.) **Stethoscope**

Put on the stethoscope allowing the earpieces of the head frame to fit well inside your ears. These arms can be rotated within their spring so that they can be positioned individually.

7.) Now place the arm with the cuff back on the table, comfortably, slightly angled, resting.

Take the manometer into your right hand, holding it around the pump, with the thumb and forefinger on the valve screw. Close the valve with a right turn and pump up the cuff steadily, watching the manometer.

Continue to pump until the indicator on the manometer exceeds your usual blood pressure by ca. 30 mmHg.

Should you still hear pulsing sounds through the stethoscope, continue to inflate until these sounds stop, which indicates that the circulation is disrupted.

8.) Open the valve very slightly with thumb and forefinger and watch the slow decrease in pressure on the manometer. The rate of deflation between systole and diastole should be in the range of 2 – 3 mmHg per second, ie there should be an interval of ca. 20 – 25 seconds for the pointer on the scale to move from 150 to 100.

- **When you hear the first clear arterial sound (knocking sound) through the stethoscope you can read the systolic (upper) blood pressure on the manometer scale.**
- **When the arterial sound stops, you can read the diastolic (lower) blood pressure.**

9.) After you have taken your readings, open the valve screw fully turning it to the left, to allow all air to be released swiftly from the cuff.

If you want to repeat your measurement, allow at least two minutes for your circulation to get back to normal. If no further measurements are required, remove the cuff.

Technical Data

Measurement range: 0 – 300 mmHg

Accuracy of Pressure shown: ± 3 mmHg

Cleaning of Unit and Cuff

To clean the unit only use a soft dry cloth.

To clean the cuff, first remove the rubber bag.

The cuff can be washed by hand at maximum 30 degrees Celsius.

Guarantee information

For details on the conditions of guarantee, please refer to the guarantee card supplied.

A claim for guarantee can only be allowed if the Certificate of Warranty, completed and stamped by the dealer, is enclosed with the instrument.

For both guarantee and repair work, please send the instrument carefully packed and postage paid either to your authorized dealer or directly to:

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BOSCH + SOHN GmbH u. Co. KG