

Discover our new products: bosolux auris and bosolux penlight

**BOSCH
+SOHN** **boso**

bosolux auris | Pocket-sized otoscope

Medical product - Made in Germany

- LED ring with 6 LEDs (concentric), dimmable
- Maintenance free – no lamp replacement necessary
- Swivel-mounted 3-piece magnifying lenses
- Insufflation bulb for pneumatic test
- Automatic shut-off after 3 minutes
- Incl. alkaline manganese batteries
- Incl. 20 single-use ear specula (10 each)
in sizes of 2.5mm + 4.0 mm

ACCESSORY PARTS

- Single-use ear specula of 2.5mm or
4.0mm in 100 piece or 1.000 piece
- Insufflation bulb for pneumatic testing
- Alkaline manganese batteries 1.5 V Micro (AAA)
(1 set of 4 pieces)



bosolux penlight | Diagnostic light for eye and throat diagnostics

Medical product - Made in Germany

- High power LED for eye diagnostics
- Housing made of high-quality ABS
- Light intensity: 25.000 – 30.000 mcd
- Incl. alkaline manganese batteries

ACCESSORY PARTS:

- Combination holder with laryngeal mirror
- Blue filter
- Alkaline manganese batteries 1.5 V Micro (AAA)
(1 set of 4 pieces)

